
**Telea**  
MEDICAL

Q-Touch

# Q-Touch

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QMR (Quantum Molecular Resonance) is an innovative technology that involves the interaction of molecular bonds that make up tissue with an electric field, having specific spectrum of frequencies transferred simultaneously [1].

In the context of Physiotherapy, the QMR-stimulation stands out from the more common radiofrequency therapies:

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

contrazioni e decontrazioni successive delle cellule [2];  
series of subsequent contractions and relaxations of cells [2];

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release of calcium ions that activates the intracellular  
pathways calcium-dependent [2].

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The combination of these effects, by inducing cells proliferation (growth of new cells), enables the regeneration of QMR-stimulated tissues [2].

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
In order to ensure maximum safety for the patient, Telea has performed a series of experiments with the purpose to ascertain that the QMR-stimulation is totally safe and it does not induce any type of cellular damage [3] [4] [5].

Scientific studies conducted on genes, chromosomes, proteins and cellular apoptosis are leaded by Laboratorio di Terapia Cellulare Avanzata (Vicenza Hospital), CIBIO and Department of Human Anatomy and Physiology (University of Padua).

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“Viability was not affected by QMR; indeed, more than 95% of cells were alive similarly to the controls, with low variability between the different MSC batches and settings. These result confirmed the morphological observations. MSC proliferation was not affected by QMR showing no significant differences between controls and QMR-treated samples at different setting and times” [3].

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“In any case no signs of apoptosis were observed. This suggests that the application of high frequency currents do not have harmful effects on the cells (myoblasts and fibroblasts)” [4].

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“The karyotype did not show the onset of chromosomal damage, such as rearrangements or breakages compared to the control” [5].

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Definitely “Phenotype and multilineage differentiation were maintained and no alteration of cellular viability and proliferation were observed between treated msc samples and control” [3].

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


Evident absence of cellular damage

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Adipogenic and osteogenic differentiation after QMR cycles of stimulation.

Panels display one representative experiment showing the final outcome in MSC multilineage differentiation after 21 days of induction.

QMR-treated (80 nominal powers) and untreated samples (CTL+) were induced to differentiation. Osteogenic differentiation (A)

after QMR stimulation was assessed using Alizarin Red. Adipogenic differentiation (B) were detected using Oil Red O.

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The perfect combination of the Therapist's expertise and the QMR-stimulation

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The innovative Q-Touch gloves enable to deliver QMR-energy straight from Therapist's hands to the patients, without, of course, involving the Operator itself.

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L'efficacia  
della terapia QMR  
QMR-stimulation

=

UNITI IN UN UNICO  
POTENTE STRUMENTO  
JOINED IN AN UNIQUE  
AND EFFECTIVE TOOL

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Only the use of QMR (patented technology – highest frequencies, lowest power) enables the energy to overcome the Q-Touch gloves. No other device available in the global market, can provide this QMR- treatment.

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The energy is delivered by the device and reaches the patient making a circuit that closes itself on a floating neutral plate.  
The involvement in this QMR-therapy exclusively occurs when the patient is in contact with the neutral floating plate.  
The operator, although wearing Q-Touch gloves, delivers QMR-energy that does not involve itself, because the therapist is outside from the electrical circuit.

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
Further and greater safety for operators, by using Q-Touch gloves, has been observed attesting the absence of leakage current.

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The graph shows the quantity of power perceived by the patient and the Operator at different display values.

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
### Comparison between power and regulation

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The measurements made show that during the QMR-therapy there is no power absorbed by the Operator at any regulation.

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Connect the neutral plate to the device and place it under the patient.

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Connect gloves cable to the specific connector.

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Wear Q-Touch gloves with the metallic connector facing up, after removing possible rings, bracelets.

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Wear up vinyl disposable gloves with high wrist. These gloves must cover up the Q-Touch gloves wristbands.

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Connect cables to the metallic connector of Q-Touch gloves.

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7

8 Impostare l'energia QMR da erogare.

Set up the QMR-energy to be delivered.

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Pour a sufficient amount of slow-absorption neutral emulsion to facilitate the sliding along treated areas.

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
Depending from the manipulation technique and, consequentially, from the dimensions of the surface in contact with the patient, the intensity of the warm sensation felt by the patient may vary accordingly.

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

Thermal photos have been taken with thermal infrared-camera (FLIR E60), during the manipulation with Q-Touch shown in normal photos.

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The operator, by using as main tool its hands, reaches and treats with QMR-therapy any part of the body, even areas where adhesive electrodes and/or handpiece find a more difficult utilization.

For example elbows, shoulders, internal parts among fingers, ankles, cervical areas etc.

Q-Touch enables precision, effectiveness and results are reached without difficulties.

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The patient will wear Q-Touch gloves over vinyl single-use gloves and, by being in contact with the neutral floating plate, she/he will be involved within the circuit of QMR-energy.

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Telea Electronic Engineering S.r.l.